


mindray

specialty transducers



CB10-4 Simultaneous Bi-Plane Endo-Cavitary Transducer

Frequency 9-4 MHz, Harmonic Frequency 10-8 MHz, Doppler Frequency 5-4.4 MHz. FOV 160 degrees. Fully supports 2D, Colorflow, Power Doppler, PW Doppler. S & T Plane views.



6LB7 Radial/Linear Bi-Plane Transducer Endo-Cavitary Transducer

Frequency 8-5 MHz, Harmonic Frequency 10-8 MHz.
Doppler Frequency 5.7-4.4 MHz.
Fully supports 2D, Colorflow, Power Doppler,
PW Doppler.



6LE7 Linear Endo-Cavitary Transducer

Frequency 8-5 MHz, Harmonic Frequency 10-8 MHz.
Doppler Frequency 5.7-4.4 MHz.
Fully supports 2D, Colorflow, Power Doppler,
PW Doppler.

for additional information, contact:

Imaging Concepts, Inc. 2818 Hungary Road Richmond, VA 23228 800-228-0060