

# M-5 Laptop Ultrasound System

The model M-5 is a new extraordinary laptop ultrasound system designed for mobile sonography. This system has phased array capability for cardiac echos, excellent resolution for general abdominal, OB/GYN and urology imaging as well as superb Small Parts and Vascular imaging. The unit comes standard with a large hard drive for image and report storage, has exquisite image and report management capability and can be configured with any standard Windows printer. The M-5 has a large assortment of transducers for most any diagnostic ultrasound test. The system promotes features such as i-Clear for excellent greyscale imaging, THI-tissue harmonic imaging, FTO-fine tissue optimization, Quadrant imaging as well as many other fine features. The system has 2D, Color-flow, Power Doppler and PW Doppler capability.



M-5 with Optional Cart



M-5 Standard Carry Case



M-5 Stand Alone

For Information, contact:

***Imaging Concepts, Inc.***

***2818 Hungary Road***

***Henrico, Virginia 23228***

***sales@imagingconcepts.com***

***800-228-0060***