

Extreme Resolution, Portable **U-5** Laptop Ultrasound System *from Medison*

New from Medison...

A portable Ultrasound System with
Extraordinary Resolution.

This product has capability for all
scanning modalities, 2D, 3D, 4D,
Colorflow, PW Doppler, Power Dop-
pler and Cardiac.

Small, light weight and powerful, the
U-5 can be used in all markets, includ-
ing OB, GYN, Urology, Small Parts,
Breast, Cardiology, Pediatrics, Radi-
ology, Interventional, ER and MSK.

This unit has comprehensive mea-
surement packages and it is DICOM
and DICOM SR compliant.



L5-12

- Linear Array Type
- Broadband Frequency Range: 5MHz - 12MHz
- **Application: Small Parts, Breast, Vascular, Musculoskeletal**

C3-7

- Curved Array Type
- Broadband Frequency Range: 3MHz - 7MHz
- **Application: Abdomen, Gynecology, OB, OB Early, Fetal Heart**

P2-4

- Phased Array Type
- Broadband Frequency Range: 2MHz - 4MHz
- **Application: Cardiology, TCD, Neonatal**

EV4-9

- Curved Array Type
- Broadband Frequency Range: 4MHz - 9MHz
- **Application: Gynecology, OB, OB Early, Fetal Heart, Urology**

3D 2-6

- Curved Array Type
- Broadband Frequency Range: 2MHz - 6MHz
- **Application: Abdomen, Gynecology, OB, OB Early, Fetal Heart**

The U-5 comes with a full range of scanning transducers to choose from for your specific application.

Also available is a convenient cart that is lightweight and very portable. The cart can also be adjusted for your specific height requirements.

The U-5 can be operated by AC power or can be disconnected and operated on the included internal battery.

This makes the U-5 an excellent product for any application where true portability is important.

Call us toll free today to see this exciting new, high resolution product, that rivals any high end system...

Medison U-5

Distributed by: Imaging Concepts, Inc.
2818 Hungary Road
Henrico, Virginia 23228

800-228-0060